

BAILEY

MAT TABLES

PROFESSIONAL HI-LOW: ELECTRIC

Mat tables have been specially designed to provide an ample exercise area at the height most convenient to wheelchair and ambulatory people. Upholstered top tables may be joined to form a larger exercise area. Mat tables are made with the same quality construction as our treatment tables and can be protected from wheelchair and prostheses damage with an optional plastic edge covering. Both upholstered tops and accessory mats for raised rim tables are 2" polyurethane foam covered with nylon reinforced Herculite 20 material. Standard color is Navy. Green, Yellow, White, and Black are available at no extra charge. Herculite 80 and Vinyl materials are also available at no extra charge. Natural finish is standard. Walnut finish is available.

The Professional Mat Table is designed with an advanced linear motor lift system. This product has all the durability and quality features that a professional demands for an efficient and economical department. A Foot Switch is standard.

- electrical height adjusts from 23" to 32" high
- 110V 60Hz
- UL hospital grade plug
- power cord 120" long
- 500 LB weight capacity



UPHOLSTERED

MODEL 9620 (4' X 7')

MODEL 9621 (5' X 7')

MODEL 9622 (6' X 8')

RAISED RIM

MODEL 9610 (4' X 7')

MODEL 9611 (5' X 7')

MODEL 9612 (6' X 8')



ACCESSORIES

Made to fit Raised Rim Mat Tables, the extra big sizes are ideal for floor exercises. Two handles on one long side make lifting and hanging very easy. Zipper closure.

MODEL 51 (4' X 7' X 2")

MODEL 53 CENTER SEAM (5' X 7' X 2")

MODEL 54 CENTER SEAM (6' X 8' X 2")

