

BAILEY

SPORTS MEDICINE

CABINET TABLES:

MODEL 4310 & 4311

- 30" tall x 30" wide x 78" long
- 2" thick vinyl upholstery with natural or dark walnut wood finishes
- twin sliding doors on hardwood tracks with enclosed storage area (MODEL 4310)
- four doors with wire fingerpulls (classic style) or optional recessed pulls on pro style doors (MODEL 4311)
- heavy duty adjustable legrests and backrest with 8 locking positions (up to 65 degrees for backrest, and 32 degrees for legrests)
- weight capacity- 350 lbs
- ships fully assembled- ship weight 260 lbs

**MODEL 4310****MODEL 4311**

CLASSROOM TRAINING TABLE:

H Brace Treatment Table: MODEL 487

- solid hardwood and birch plywood construction
- natural finish
- 22" chair clearance standard (optional heights available)
- "English Oak" laminated top with radius corners
- 2" thick vinyl covered mat with high-density foam (removable with hook & loop attachment)
- table size 30" tall x 30" wide x 78" long x 30"H
- requires on site assembly

**MODEL 487**

PROFESSIONAL EXERCISE TABLE:

Quadriceps Knee Tables: MODEL 345 & 350

- 30" tall x 30" wide x 22" deep
- 1" thick vinyl upholstery
- twin torque units use standard 1" ID disk weight plates (not included)
- 0-90 degree adjustable back support and headrest (MODEL 345).
- twin foam padded handles for ease of use (MODEL 345)
- adjustable thigh strap for stabilization (MODEL 345 or MODEL 3901 on MODEL 350)
- twin adjustable weight slides & stationary weight posts
- twin adjustable slide shin bolsters
- full range locking torque dials
- torque unit capacity 37.5 lbs, weight capacity 350 lbs
- adjustable "click lock" strut

**MODEL 345****MODEL 350**

TREATMENT TABLES:

**MODEL 441****MODEL 432**

WHIRLPOOL TABLES:

**MODEL 7724****MODEL 7720**

BAILEY MANUFACTURING CO.

P.O. BOX 130 LODI, OHIO 44254-1030

IN USA: PH 800-321-8372 FAX: 800-224-5390

OUTSIDE USA: PH 330-948-1080 FAX 330-948-4439

WWW.BAILEYMFG.COM

bailey_mfg_lodi



Bailey Manufacturing Company

CONTACT: