

# BAILEY

## SPORTS MEDICINE

### TAPING TABLES

*Taping Tables are designed to make the taping of athletes quick and efficient. Built-in shelves and cabinets put tape, scissors and other supplies within easy reach. And the compact size is just right for smaller rooms.*

#### Single Taping Station: MODEL 27

- height choices- 36" or 42"
- 1" thick vinyl upholstery
- 32" deep x 26" wide taping tops w/ standard 28 quart wastebin supply top
- standard dual adjustable interior shelves and triple drawers
- optional locks for doors and drawers
- base size 66" wide x 30" deep
- finished all-around, use as an island or inline
- weight capacity 500 lb.
- ships fully assembled

\*Model 27 may be daisy chained, with or without Model 26, to make any combination of space saving taping stations!

#### Dual Taping Station: MODEL 26

- height choices- 36" or 42"
- 1" thick vinyl upholstery
- 32" deep x 26" wide taping tops w/ standard 28 quart wastebin supply top
- full depth adjustable interior shelves and triple drawers
- optional locks for doors and drawers
- base size 66" wide x 30" deep
- finished all-around, use as an island or inline
- weight capacity 500 lb. per seat

**MODEL 26****MODEL 27**

#### Deluxe Taping Table: MODEL 25

- height choices- 30", 36" or 42"
  - 1" thick vinyl upholstery
  - 32" deep x 26" wide taping tops w/ standard 28 quart wastebin supply top
  - natural or dark walnut wood finishes
  - full rear shelf and end storage shelves
  - stationary 60" back support
  - weight capacity 400 lb.
- (shown with custom embroidery option)

#### Fully Enclosed Taping Table: MODEL 10

- easy to clean 1 1/4" wood grain laminate top
- 3/4" multi-layered birch sides and shelves
- standard height- 36"
- optional height of 30" or 42"
- 26" wide x 40" long
- natural finish is standard, walnut finish is available
- upholstered top available as a special request item



#### End Shelf Taping Table: MODEL 12

- standard height- 36"
- optional height of 30" or 42"
- 26" wide x 40" long upholstered top
- stationery end shelves for storage
- natural finish is standard, walnut finish is available
- laminated tops are available as a custom option

#### Taping Table: MODEL 14

- standard height- 36"
- optional height of 30" or 42"
- 26" wide x 48" long upholstered top
- natural finish is standard, walnut finish is available



### BAILEY MANUFACTURING CO.

P.O. BOX 130 LODI, OHIO 44254-1030

IN USA: PH 800-321-8372 FAX: 800-224-5390

OUTSIDE USA: PH 330-948-1080 FAX 330-948-4439

[WWW.BAILEYMFG.COM](http://WWW.BAILEYMFG.COM)



bailey\_mfg\_lodi



Bailey Manufacturing Company

### CONTACT:

*Over 60 years of experience in crafting equipment for physical therapy, occupational therapy, work hardening, and sports medicine professionals. Our products are sold through a national network of authorized dealers. Contact us for the dealers nearest to you.*